

First Aid

Features

- Emergency eye irrigation for any work environment
- Sterile 0.9% saline solution ideal for irrigating eyes
- Suitable for workplaces with a risk of eye contamination from debris or chemicals
- Single use bottles for use as a substitute for clean water
- Size is 500ml

RS PRO HypaClens Emergency Eye Wash Bottle, 500ml

RS Stock No.: 787-3521



RS PRO is the own brand of RS. The RS PRO Seal of Approval is your assurance of professional quality, a guarantee that every part is rigorously tested, inspected, and audited against demanding standards. Making RS PRO the Smart Choice for our customers.

First Aid

Product Description

The HypaClens Sterile Eye Wash is ideal for flushing eyes that have been contaminated with particles, smoke, dust or liquid. The bottles contain a sterile 0.9% saline solution that gently cleanses the eye.

HSE guidelines state that at least 1 litre of eye wash should be available at any given time in work environments where flying debris or chemicals could pose a danger to the eye. Following eye contamination with chemicals or small particles, you should wash the eye out with either clean water or sterile fluid from a sealed container for at least 10 minutes.

The HypaClens Sterile Eye Wash is intended as a substitute for clean running water in first aid situations – each bottle should only be used once. After the initial use of the bottle, its contents can no longer be considered sterile and could pose a risk of infection. After use, any excess fluid remaining in the bottle should be discarded.

The HypaClens Sterile Eye Wash is also available in 250ml.

General Specifications

Product Dimensions	19x24x7.5 cm
Product Weight (grams)	545
Wall Bracket Included	No
Medical Device Class	Class I Sterile
Application	First Aid
Suitable for	rinsing foreign particles out of eyes
Shelf Life	Up to 36 months

Classification

HS Code	3004900000
---------	------------